



Active Policy Name: Basic ALC Profile ID: basic Creator: DRG

- SOP UID
- Other UID
- Study UID
- Series UID
- Remove Institution Data
- Anonymise Patient Data
- Patient ID
- Remove Patient Data
- Ambiguate Date of Birth
- Remove Other

Transformer: FieldDestroyer Select

Parameters:



DICOM Tags

Targets

- (0000,0000) GenericGroupLength UL 1
- (0000,0002) AffectedSOPClassUID UI 1
- (0000,0003) RequestedSOPClassUID UI 1
- (0000,0100) CommandField US 1
- (0000,0110) MessageID US 1
- (0000,0120) MessageIDBeingRespondedTo US 1
- (0000,0600) MoveDestination AE 1
- (0000,0700) Priority US 1
- (0000,0800) DataSetType US 1

- PatientsBirthTime
- PatientsSex
- OtherPatientIDs
- PatientsAge
- PatientsSize
- PatientsWeight
- PatientsAddress
- PatientsMothersBirthName
- MedicalRecordLocator
- PatientsTelephoneNumbers

